

continued from page

Testing Spectroscopy Files → 5/16/2024

Run 1 → Intensity Measurement with IPA @ Room Temp Deuterium

Run 2 → Transmission View, IPA @ Room Temp, Deuterium

↳ Data Update Rate 5 ms; Integration Time 15 ms; ~~Scans~~

Scans to Average 5 Boxcar width 5

Kevan Williams William Mulberry

Run 3 → Intensity with IPA @ ~~Room Temp~~ Deuterium @ 42°C

Run 4 → Intensity with IPA Deuterium @ 65°C

Run 5 → Intensity with IPA Deuterium @ 74°C

-Started boiling in this range.

Kevan Williams